

## HOW TO PUT ON YOUR NEW PROSTHESIS



Roll liner completely inside out before rolling on. Make sure longer side is on outside.



Make sure it's rolled evenly all the way up with no wrinkles.



Apply socks as needed to tighten the fit.



Make sure to thread the KISS strap through the bottom hole. Then pull them all the way up.



Locate the slot inside of your prosthesis and feed the strap through, Velcro side down.



Push the strap through until there is enough to grab and pull the rest out.



Pull the strap up as tight as possible and thread it through the metal loops.



Pull the strap down tightly and fasten the Velcro.



Once standing the strap can be further tightened to provide a more secure fit.

**FOR QUESTIONS OR CONCERNS PLEASE CALL OUR OFFICE (732) 380-0366**

**Manfredi Orthotic & Prosthetic Affiliates  
749 B Hope Rd, Suite C Eatontown, NJ 07724**