

How to Don a magnetic locking suspension system



Roll on Gel liner with maglock fitting, make sure you leave no air space at end



Add as many socks as necessary to make socket snug

Engage mag liner to lock and push in button and slide lock to lock position to secure



Pull up Suspension sleeve, get out all the wrinkles, and make sure the top 2 "

Is in contact with your skin

If you have any questions or concerns please call the office at **732-380-0366**

Manfredi Orthotic & Prosthetic Affiliates

749 Hope rd, Suite C Eatontown, NJ 07724

manfredioandp.com