

# How to put on a prosthesis with a locking liner



Hold gel liner against your limb



Make sure there are no wrinkles



roll liner up completely



Add socks as needed making sure pin is through the hole



Engage the pin in the lock



Push down on limb until you hear a "click"



Roll up suspension sleeve



Make sure to smooth out wrinkles



To take off, push button

If you have any questions or concerns please call the office at 732-380-0366

**Manfredi Orthotic & Prosthetic Affiliates**

749 Hope rd, Suite C Eatontown, NJ 07724